

# ABDOMEN PHANTOM

Age  
Category


Adult

Body  
Region

Abdomen

Target  
Modality

CT



This abdomen phantom can be used in CT for image quality evaluation and for training.

The phantom simulates a contrast medium enhanced abdomen in late arterial phase.

The phantom provides a detailed and realistic simulation of soft and bone tissue, including small details such as lymph nodes. Air voids are filled with a cellulose-polymer composite of approx. -80 HU.

The phantom can be used for common methods of image quality evaluation such as visual grading analysis or contrast-to-noise measurement.



# ABDOMEN PHANTOM



## Specifications

Size	Approx. 267 x 178 x 167 mm
Weight	Approx. 5700 g
Base material	Cellulose-polymer composite
Optimal tube voltage	120 kVp (cf page 3) - adaptable upon request -

---

## Similar products

- Abdomen phantoms with liver lesions
- Abdomen phantoms with pancreatic lesions

For more information visit  
[www.phantomx.de](http://www.phantomx.de)

# ABDOMEN PHANTOM

## General indications


- The phantom is made of a cellulose-polymer composite material with properties similar to hardwood. If handled carefully, it will last a long time.
- The phantom is coated with a protective layer. If the protective layer is undamaged, the phantom can be cleaned using a damp cloth (water or mild detergent).
- Protect from direct sunlight.
- Maintain a storage temperature of 10 °C to 30 °C. If the phantom is exposed to temperatures below -10 °C or above 45 °C, it can be severely damaged.
- The phantom is not equipped for dose measurements with dosimeters and it is not suited for material characterization with dual energy CT.
- The phantom is not certified as medical device.
- Air voids are filled with cellulose-polymer composite of approx. -80 HU.

Phantom based on modified data, originally from "Abdominal\_skin\_test01 1.0.0" by EranJ licensed under CC BY 4.0.


## Attenuation properties

### Soft Tissue

Linear attenuation coefficients [ $\text{cm}^{-1}$ ] (calculated)




### Hounsfield units (calculated)




### Bone Tissue

Linear attenuation coefficients [ $\text{cm}^{-1}$ ] (calculated)



### Hounsfield units (calculated)



Tissue Reference: Woodard HQ, White DR. The composition of body tissues. Br J Radiol. 1986.